

# version 24.20.0 Release Notes



## WHAT'S NEW

### Zoetis Reference Labs

#### Add-ons to Tests

Operators can add onto tests that have already been submitted to Zoetis. The tests will have to be mapped and configured (using the product attachments and the lab request sheets) to generate lab requests for Zoetis for the add-ons to work.

#### Add-ons from ZoetisDx

Operators can place an order from Impromed, view the results within ZoetisDx, then add additional tests onto the order. When they do this, updates are sent to Impromed on what was added and allowing the operator to reconcile the additional tests in Impromed so they can still do the charge capture.

# WHAT'S NEW

## All Lab Integrations

### Defer Submitting Lab Requests

Operators will be able to choose “Submit Later” when they are invoicing lab items which trigger lab requests. Operators can choose to submit when they continue to add products which trigger the generation of lab requests or when the invoice is concluded. *Note: This feature is only available when using the product attachments to trigger lab requests when adding items to the invoice.*

### Orphan Results Charge Capture

All lab partners will have the ability to capture charges for orphaned lab results during the import process.

# ADDITIONAL RELEASE ITEMS

## Covetrus Comms

413289	The option to use Covetrus Comms for email will not be checked by default.
--------	--

## Dashboard

408244	The Reminder Status Report has been corrected to display the word “Performed”.
410953	The list of Doctor Preference will match between the Dashboard and the Database module.

## Lab Integrations

410586	The providers listed in the Lab Integrations filter will show only those employees that can add or modify labs.
413037	Corrected the legacy unidirectional lab travel sheet so it will not display the error that the product is not available when the product is available.
415267	Zoetis lab requests creating intermittent errors has been corrected.
415268	The “Show Web Order” will be hidden for any lab request with a status marked “Not Complete” that has a temporary ID.
416417	Impromed allows labs from different partners to be imported if their Lab IDs are identical to AXIS-Q.
412709	When users add products that generate a lab request for more than one patient on the same invoice (before clicking OK), all the requisitions made will display correctly.

# ADDITIONAL RELEASE ITEMS

## Medical Records

---

<b>90147</b>	Corrected the error when transferring a patient or medical records.
--------------	---

---

## Mobile

---

<b>415843</b>	The Mobile Sync will recognize when there is one user enabled and ignore any disabled users with the same name so it does not cause an error.
---------------	---

---

## Reports

---

<b>414159</b>	The NE 5.0 Control Substance log export has been updated to include all required fields during submission.
---------------	--

---