

SETTINGS FOR DYMO PRINTERS

Settings For DYMO® Printers

1. In AVImark, go to **Utilities | Printer Setup**.
2. Select the **DYMO** for Drug Labels and select the **Quality** as Label Printer.
3. Click the **Form Assignment** tab and select **30258 Diskette** in the box to the right of the DYMO, if you are using the recommended label size of 2 1/8 x 2 3/4.
4. Click **Done** and close AVImark (*don't minimize the program*).
5. Depending upon your version of Windows, open the **Control Panel**. Search for **Printers** then view **Devices and Printers**.
6. Right-click on the **DYMO** printer and click **Printer Properties**.
7. Click **Preferences**.
8. For **Orientation**, select the **Landscape** option.
9. Click **Advanced** located on the right.
10. Select **30258 Diskette** as the paper size.
11. Click **OK**,

Follow the steps listed here on every computer that has a DYMO printer installed on.